

SPECT Reconstructed Spatial Resolution without Scatter



SPECT Reconstructed Spatial Resolution without Scatter (nmtk-rsr)

This fixture follows NEMA Nu-1 System Alignment and SPECT Reconstructed Spatial Resolution without Scatter test equipment requirements.

Items included:

- Aluminum base plate (20cm x 46cm) – slotted for three capillary glass tubes (75mm separation with 50mm axial offset) Supply of 20 cm long (1.0 mm ID) glass capillary tubes and sealer
- Supply of 20 cm long 1.0mm diameter (ID) glass capillary tubes and sealer

QA Benchmark, LLC
1101 Resort Road #113 | Ellicott City, MD 21042
p. 443.980.6349
e: qabenchmark2017@gmail.com
w: qabenchmark.com

